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Dr. Rostislav Boltyanskiy,
Spheryx Senior Scientist II & Director of Scientific Communications

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info@shoko-sc.co.jp

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1. Krishnan M, Nalaskowski, JW, and Cook, LM. Chemical Mechanical Planarization: Slurry Chemistry, Materials, and Mechanisms. Chem. Rev. 2010, 110, 1, 178–204. doi.org/10.1021/cr900170z.
2. Brahma N, Talbot JB. Effects of CMP slurry additives on the agglomeration of alumina nanoparticles 1: General aggregation rate behavior. Journal of Colloid and Interface Science. 2014, 419, 56-60. doi.org/10.1016/j.jcis.2013.12.029.
3. Beyer, K. Guthrie WL, Makarewicz SR et al. Chem-Mech Polishing Method for Producing Coplanar metal/insulator films on a substrate. 1990. US Patent number: 4,944,836.
4. Transparency Market Research. CMP Slurry Market. 2022.
<https://www.transparencymarketresearch.com/cmp-slurry-market.html>
5. Khanna AJ, Gupta S, Kumar P, Chang FC, and Singh RK. Quantification of shear induced agglomeration in chemical mechanical polishing slurries under different chemical environments. Microelectron. Eng. 2019, 210, 1. doi.org/10.1016/j.mee.2019.03.012
6. Khanna AJ, Gupta S, Kumar P, Chang FC, and Singh RK. Study of agglomeration behavior of chemical mechanical polishing slurry under controlled shear environments. ECS J. Solid State Sci. Technol. 2018, 7, P238. doi.org/10.1149/2.0091805jss
7. Raghavan SR and Khan SA. Shear-induced microstructural changes in flocculated suspensions of fumed silica. J. Rheol. 1995, 39, 1311. doi.org/10.1122/1.550638
8. Bakier M, Suzuko K, Khojornrungruang, P. Study on Nanoparticle Agglomeration During Chemical Mechanical Polishing (CMP) Performance. Journal of Nanofluids, 2021, 10, 3, 420-430. doi.org/10.1166/jon.2021.1791.
9. Cheong FC, Kasimbeg P, Ruffner DB, Hliang EH, Blusewicz JM, Philips LA, and Grier DG. Holographic characterization of colloidal particles in turbid media. Applied Physics Letters, 2017, 111, 153702. <https://doi.org/10.1063/1.4999101>
10. Cheong FC, Kasimbeg P, Winters A, Blusewicz JM, Hliang E, Ruffner DB, Ihnfeldt RV, Grier DG, Philips LA. Holographic characterization of agglomerates in CMP slurries. ICPT 2018, 216.
11. Kosier S, Blackford D. A New Method for Determining the Size Distribution of Particles in CMP Slurries. ICPT 2018, 120.
12. Kiernan J, Bengochea LV, Johnston R, Kavaljer M. In-Line Refractive Index in Assay Characterization of Incoming Fresh and Effluent Spent CMP Slurry, ICPT 2018, 132.

13. Al-Gebory L, Menguc AP. The effect of pH on particle agglomeration and optical properties of nanoparticle suspensions. *Journal of Quantitative Spectroscopy and Radiative Transfer*. 2018, 219, 46-60. doi.org/10.1016/j.jqsrt.2018.07.020
14. Godymchuk A, Papina I, Karepina E, et al. Agglomeration of iron oxide nanoparticles: pH effect is stronger than amino acid acidity. *Journal of Nanoparticle Research*. 2019, 21, 208. doi.org/10.1007/s11051-019-4634-y
15. Loosli F, Stoll S. Effect of surfactants, pH and water hardness on the surface properties and agglomeration behavior of engineered TiO₂ nanoparticles. *Environmental Science: Nano*. 2017, 1. doi.org/10.1039/C6EN00339G
16. Sadiq SA, Waseem Z, Hanif S et al. Investigation of the role of pH on structural and morphological properties of titanium dioxide nanoparticles. *OP Conf. Ser.: Mater. Sci. Eng*. 2020, 863. doi.org/10.1088/1757-899X/863/1/012046
17. Palla BJ, Shah DO. Stabilization of High Ionic Strength Slurries Using the Synergistic Effects of a Mixed Surfactant System. *J Colloid Interface Sci*. 2000, 223, 1, 102-111. doi: 10.1006/jcis.1999.6665
18. Lagudu URK, Isono S, Krishnan S, Babu SV. Role of ionic strength in chemical mechanical polishing of silicon carbide using silica slurries. *Colloids and Surfaces A: Physicochemical and Engineering Aspects*. 2014, 445, 119-127. doi.org/10.1016/j.colsurfa.2014.01.038
19. Chang FC, Tanawade S, Singh RK. Effects of Stress-Induced Particle Agglomeration on Defectivity during CMP of Low-*k* Dielectrics. *Journal of Electrochemical Society*. 2008, 156, H39. doi.org/10.1149/1.3005778
20. Khanna AJ, Chang FC, Gupta S, Kumar P. Characterization of the nature of shear-induced agglomerates as hard and soft in chemical mechanical polishing slurries. *Journal of Vacuum Science & Technology B*. 2019, 011207. doi.org/10.1116/1.5065516
21. Chang FC, Kumar P, Singh R, et al. Role of interparticle forces during stress-induced agglomeration of CMP slurries. *Colloids and Surfaces A: Physicochemical and Engineering Aspects*. 2014, 445. doi.org/10.1016/j.colsurfa.2011.09.001