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1. Vishtal, A. and A. Kraslawski, Challenges in industrial applications of technical lignins. *Bioresources*, 2011. 6(3): p. 3547-3568.
2. Lauten, R.A., B.O. Myrvold, and S.A. Gundersen, New Developments in the Commercial Utilization of Lignosulfonates, in *Surfactants from Renewable Resources* 2010, John Wiley & Sons, Ltd. p. 269-283.
3. Lora, J., Chapter 10 - Industrial Commercial Lignins: Sources, Properties and Applications, in *Monomers, Polymers and Composites from Renewable Resources*, M.N.B. Gandini, Editor 2008, Elsevier: Amsterdam. p. 225-241.
4. Marques, A.P., et al., Chemical Composition of Spent Liquors from Acidic Magnesium-Based Sulphite Pulping of Eucalyptus globulus. *J. Wood Chem. and Tech.*, 2009. 29(4): p. 322-336.
5. Mika, M. and N. Marianne, Utilization of Membrane

Processes in Treating Various Effluents Generated in Pulp and Paper Industry, in *Handbook of Membrane Separations* 2008, CRC Press. p. 981-1006.

6. Albinsson, B., et al., The origin of lignin fluorescence. *J. Mol. Struct.*, 1999. 508(1-3): p. 19- 27.
7. Fredheim, G.E., S.M. Braaten, and B.E. Christensen, Molecular weight determination of lignosulfonates by size-exclusion chromatography and multi-angle laser light scattering. *J. Chromatography. A*, 2002. 942(1-2): p. 191-199.
8. Sulaeva, I., et al., Molar Mass Characterization of Crude Lignosulfonates by Asymmetric Flow Field-Flow Fractionation. *ACS Sustainable Chem. Eng.*, DOI: 10.1021/acssuschemeng.8b02856
9. Ringena, O., et al., Size-exclusion chromatography of technical lignins in dimethyl sulfoxide/water and dimethylacetamide. *J. Chromatography A*, 2006. 1102(1-2): p. 154-163.