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Dr. Rostislav Boltyanskiy,
Spheryx Senior Scientist II & Director of Scientific Communications

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1. Schellekens H. Bioequivalence and the immunogenicity of biopharmaceuticals. *Nat Rev Drug Discov.* 2002;1:457– 462.
2. Wang W. Protein aggregation and its inhibition in bio- pharmaceuticals. *Int J Pharma.* 2005;289:1–30.
3. Carpenter JF, Randolph TW, Jiskoot W, Commelin DJA, Middaugh CR, Winter G, et al. Overlooking sub- visible particles in therapeutic protein products: Gaps that may compromise product quality. *J Pharm Sci.* 2009;4:1201– 1205.

4. Singh SK, Afonina N, Awwad M, Bechtold-Peters K, Blue JT, Chou D, et al. An industry perspective on the monitoring of subvisible particles as a quality attribute for protein therapeutics. *J Pharm Sci.* 2010;99:3302–3321.
5. den Englesman J, Gandel P, Smulders R, Koll H, Smith B, Bassarab S, et al. Strategies for the assessment of protein aggregates in pharmaceutical biotech product development. *Pharm Res.* 2011;28:920–933.
6. Ratanji KD, Derrick JP, Dearman RJ, Kimber I. Immunogenicity of therapeutic proteins: influence of aggregation. *J. Immunotoxicol.* 2014;11:99-109.
7. Food and Drug Administration. Guidance for Industry: Immunogenicity Assessment for Therapeutic Protein Products. U.S. Department of Health and Human Services; 2014.
8. Panchal J, Kotarek J, Marszal E, Topp EM. Analyzing subvisible particles in protein drug products: a comparison of Dynamic Light Scattering (DLS) and Resonant Mass Measurement (RMM). *AAPS J.* 2014;16(3):440– 451.
9. Hamrang Z, Hussain M, Tingey K, Tracka M, Casas-Finet J, Uddin S, et al. Characterization of stress induced aggregate size distributions and morphological changes of a bi-specific antibody using orthogonal techniques. *J Pharm Sci.* 2015;104(8):2473–2481.
10. Sharma DK, King D, Oma P, Merchant C. Micro-Flow Imaging: Flow microscopy applied to subvisible particulate analysis in protein formulations. *AAPS J.* 2015;12(3):455–464.
11. Zoëlls S, Wiggenghorn M, Winter G, Friess W, Jiskoot W. Flow imaging microscopy for protein particle analysis – A comparative evaluation of four different analytical instruments. *AAPS J.* 2013;15(4):1200–1211.
12. Felsovalyi F, Janvier S, Jouffray S, Soukiassian H, Mangiagalli P. Silicone-oil- based subvisible particles: their detection, interactions, and regulation in prefilled container closure systems for biopharmaceuticals. *J Pharm Sci.* 2012;101:4569-4583.
13. Shah M, Rattray Z, Day K, et al. Evaluation of aggregate and silicone-oil counts in pre-filled siliconized syringes: an orthogonal study characterising the entire subvisible size range. *Int J Pharm.* 2017;519:58-66.
14. Martos A, Koch W, Jiskoot W, et al. Trends on analytical characterization of polysorbates and their degradation products in biopharmaceutical formulations. *J Pharm Sci.* 2017;106:1722-1735.
15. Scherer TM, Leung S, Owyang L, Shire, SJ. Issues and challenges of subvisible and submicron particulate analysis in protein solutions. *AAPS J.* 2012;14:236-243.
16. Wang C, Zhong X, Ruffner DB, et al. Holographic characterization of protein aggregates. *J Pharm Sci.* 2016;105:1074-1085.
17. Kasimbeg PN, Cheong FC, Ruffner DB, Blusewicz JM, Philips LA. Holographic characterization of protein aggregates in the presence of silicone oil and surfactants. *J Pharm Sci.* 2019;108:155-161.
18. Winters A, Cheong FC, Odete MA, et al. Quantitative Differentiation of Protein Aggregates from Other Subvisible Particles in Viscous Mixtures Through Holographic Characterization. *J. Pharm Sci.* 2020; 109:2405-2412.
19. Randolph TW, Schiltz E, Sederstrom D, et al. Do not drop: mechanical shock in vials causes cavitation, protein aggregation, and particle formation. *J Pharm Sci.* 2015;104:602-611.
20. Tobias W, Volkin DB, Mahler HC. Effect of solution properties on the counting and sizing of subvisible particle standards as measured by light obscuration and digital imaging methods. *Eur J Pharm Sci.* 2014;53:95-108.