

Advion is a leader in mass spectrometry & synthesis solutions. The **expression** CMS is a high performance, compact and affordable single quadrupole mass spectrometer. Its compact size allows it to fit in space-limited laboratories for direct access and immediate results for chemists requiring mass confirmation, reaction monitoring, quality control and purity analysis.

# Real Time Reaction Monitoring of a Solution Phase Peptide Synthesis Using TLC/Compact Mass Spectrometry

## Introduction


Real-time reaction monitoring based on the **expression** compact mass spectrometer (CMS) can support the synthetic chemist to overcome synthesis challenges and optimize the reaction time on-line. Peptides of pharmaceutical interest can be readily synthesized following a rapid, continuous solution-phase synthesis strategy based on Fmoc protected amino acid building blocks<sup>1</sup>. A simple model for such a reaction is the growing of analogues of the acyl carrier protein (ACP), a component of the fatty acid synthesis pathway. Thin layer chromatography (TLC) is a simple and easy way to prepare complex reaction samples for CMS analysis following a TLC/CMS approach.

## Reaction Procedure

Fmoc-Asn-(Trt)-OH (1.0 mmol, 600 mg) and O-(Benzotriazol-1-yl)-*N,N,N'*-tetramethyluronium hexafluorophosphate (HBTU, 1.0 mmol, 400 mg) were dissolved in 10 mL dichloromethane (DCM) in a 50 mL round bottom flask and stirred. The amino acid building block was activated by adding *N,N*-Diisopropylethylamine (DIEA, 2.65 mmol, 465  $\mu$ L) and stirred for 5 min. The Glycine *tert*-butyl ester hydrochloride (H-Gly-OtBu, 0.65 mmol, 109 mg) was added and the reaction commenced for 1 h at room temperature.

The Fmoc protected dipeptide can be unblocked using tris(2-aminoethyl)amine.

## Reaction Pathway



The reaction can be continued for more additions of amino acid building blocks (not shown in the reaction pathway), before the final de-protection of the product with trifluoroacetic acid (TFA).

## Development of TLC plates

TLC plates (Merck EMD, 5534-3 Silica Gel 60 F254, 5x20 cm, 0.2 mm) were developed using a mixture of Dichloromethane/Methanol/Acetic acid (9:1:0.1). TLC plates were monitored by a hand held UV light excitation at 254 nm.


## TLC spot extraction

Developed TLC plates were further processed with a semi-automated TLC spot extractor 'TLC/MS Interface' (Camag Scientific, Wilmington, NC) for TLC/CMS analysis. The TLC card was moved underneath the interface and the target area illuminated with a laser cross. The extraction head had an oval extraction shape of ca. 3x5 mm and an extraction area of ca. 15 mm<sup>2</sup>. The head was pressed against the TLC plate using 25 psi pressure and a flow of 200  $\mu$ L/min 80/20 Acetonitrile/water 0.1 vol% formic acid was flushed through the extraction head and directly to the ESI/MS source.


## MS method settings

The **expression** Compact Mass Spectrometer (CMS) (Advion Inc., Ithaca, NY) was set to scan from *m/z* 100 to *m/z* 1200 for the 3 min acquisition time in the TLC/CMS experiments.

## Results



**Figure 1:** TLC separation after clean up of the reaction solvent showing both the Fmoc blocked (left) as well as un-blocked (right) dipeptide product from reaction Step 1 (A). TLC/CMS analysis can confirm both product forms via their respective  $(M+Na)^+$  ions at  $m/z$  732.21 (B) and  $m/z$  510.15 (C).



**Figure 2:** Monitoring the first reaction step every 10 min by placing 1  $\mu$ L of the reaction solvent on a silica TLC plate and developing the plate at the end of the reaction time. UV light showed a multitude of compounds being separated on the TLC plate (A). Subsequent TLC/CMS analysis of the spots along the 10 min lane (data not shown) showed that spot 4 represented the desired product (red arrows). The Fmoc blocked product was detected at  $m/z$  732.20, the  $(M+Na)^+$  signal. Signal intensity of the respective spots show a time curve consistent with a completed reaction at 20 min.

## Summary

- Silica gel TLC/CMS was able to identify educts, products and side products of a chemical peptide synthesis reaction
- Both Fmoc blocked as well as unblocked peptides could be detected as the  $(M+Na)^+$  ion in the mass spectrometer
- TLC/CMS could monitor the chemical reaction and showed that the reaction was complete in 20 min rather than the expected 60 min.

## Literature/Acknowledgement

<sup>1</sup>Carpino LA, Ghassemi S, Ionescu D, Ismail M, Sadat-Aalae D, Turan GA, Mansour EME, Siwruk GA, Eynon JS and Morgan B: Rapid, continuous solution-phase peptide synthesis: Application to peptides of pharmaceutical interest. *Organic Process Research and Development* 2003, 1(7), 28-37

We would like to thank Shahnaz Ghassemi from Synpure LLC for her technical assistance in the peptide synthesis as well as Camag Scientific Inc., Wilmington, NC for the generous loan of a TLC/MS interface.